
Manufactured for:
Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668 USA

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Amara Full Face Mask

Intended Use

The Amara and Amara Gel Full Face Masks are intended to provide an interface for application of CPAP or bi-level therapy to patients. The masks are for single patient use in the home or multi-patient use in the hospital/institutional environment. The masks are to be used on patients (>66lbs/30kg) for whom CPAP or bi-level therapy has been prescribed.

Note: An exhalation port is built into the mask so a separate exhalation port is not required.

Note: This mask is not made with natural rubber latex or DEHP.
Symbols

Warning or
Caution

Note

Tip

Not Made With Natural
Rubber Latex

System One Resistance
Control

Warnings:

- This mask is not suitable for providing life support ventilation.
- Hand wash prior to use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or try to seal the exhalation port. Explanation of the Warning: CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be re breathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. Explanation of the Warning: When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- Consult a physician if any of the following symptoms occur: Unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, blurred vision, drying of the eyes, eye pain or eye infections.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.

- This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Consult a physician or dentist if you encounter tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- A minimum of 3cm H₂O (hPa) must be maintained when using this mask.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.
- Do not block or seal off the entrainment valve.

Contraindications

This mask may not be suitable for use on patients with the following conditions:

recent eye surgery or dry eyes; hiatal hernia, excessive reflux, impaired cough reflex, impaired cardiac sphincter function; or on patients unable to remove the mask by themselves.

Before Use Read and Understand the

Instructions Completely


- Hand wash the mask.
- Wash your face. Do not use moisturizer/lotion on your hands or face.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Verify therapy device pressure(s).

Verify the Entrainment Valve

The entrainment valve consists of an air inlet and a flapper.


With the air flow turned off, verify that the elbow with the flapper is lying flat (Figure 2A) so that room air can flow in and out through the air inlet. Next, with the airflow on, the flapper should now cover the air inlet and air from the CPAP or bi-level device should flow into the mask (Figure 2B). If the flapper does not close or does not function properly, replace the mask.

⚠ Warning: Do not block or seal the entrainment valve.


▶ Figure 1

- A Top Headgear Straps
- B Integrated Cushion/Faceplate
- C Bottom Headgear Straps
- D Entrainment Valve Air Inlet (Do Not Block)
- E Push Button Adjuster
- F Forehead Pad Support
- G Frame with Elbow
- H Exhalation Port (Do Not Block)
- I Headgear Clips
- J Quick Release Swivel

2A

2B

▶ Figure 2

Cleaning Instructions

Hand wash the mask before first use and daily. The headgear should be hand washed weekly or as needed. The headgear does not need to be removed for daily cleaning of the mask.

1. Hand wash mask in warm water with liquid dishwashing detergent.

2. Rinse thoroughly. Air dry completely before use. Make sure the mask is dry before use. Lay the headgear flat or line dry. Do not place the headgear into the dryer.

⚠ Caution: Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

⚠ Caution: Any deviation from these instructions may impact the performance of the product.

⚠ Caution: Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.

Institutional Disinfection

For **multi-patient use** in the hospital/institutional environment, use the **Disinfection Guide** to reprocess the mask between patients. These instructions can be obtained by visiting us online at www.philips.com/respironics or by contacting Philips Respironics Customer Service at 1-800-345-6443 (USA or Canada) or at 1-724-387-4000.

Dishwashing Instructions

In addition to hand washing, the mask may be cleaned in the dishwasher once a week.

⚠ Caution: Use a mild liquid dish washing detergent only to wash the mask.

1. Remove the headgear. Do not wash the headgear in the dishwasher.

2. Separate the cushion and the faceplate and wash in the top shelf of the dishwasher.

3. Air dry. Make sure the mask is dry before use.

4. Inspect all parts for damage or wear; replace any parts that have visibly deteriorated (cracking, crazing, tears, etc). Replace the mask if the cushion becomes hardened.

Philips Respironics System One Resistance Control

Your mask when combined with a Philips Respironics System One device, provides optimal resistance compensation. This mask value is  and should be set by your provider.

Notes:

- Compare the mask to the device. See your device manual if the values do not match.
- System One is not compatible with masks that use a separate exhalation device.

Achieving the Right Fit

Before putting on the mask:

1. **Headgear clips** - Grip and twist away from the frame to disconnect.
2. **Headgear** - Loosen to a large setting.

Adjusting the mask:

1. Pull the mask over your head. Hold the mask against your face.
2. **Headgear clips** - Press into place to connect. (Figure 3).
3. **Headgear** - Pull back the tabs to evenly adjust the top and bottom straps. (Figure 4).
4. **Forehead adjuster button** - Press and hold to slide the frame toward or away from your face to adjust pressure (Figure 5). This will decrease or increase the pressure on the bridge of your nose and chin.

The mask should rest comfortably against your face. The forehead pad will touch your forehead lightly. The bottom of the mask cushion should rest just above your chin with your mouth slightly open.

Using the mask:


1. Connect the flexible tubing (included with your CPAP or bi-level device) to the elbow swivel on the mask.
2. Lie down. Turn the therapy device on. Breathe normally. Keep your mouth slightly open.
3. Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments. Some air leaking is normal.

✳ Comfort Tips


- The most common mistake is over-tightening the headgear. The headgear should fit loose and comfortable. If your skin bulges around the mask or if you see red marks on your face, loosen the headgear.
- Adjust the bottom headgear straps to reduce leaks at the sides of your mouth.

Removing the Mask


Disconnect the headgear clips and pull over the top of your head (Figure 6).


▶ Figure 3

▶ Figure 4

▶ Figure 5


▶ Figure 6

Disassembly of the Mask

1. **Headgear** - Remove from the forehead pad support and headgear clips.
2. **Forehead pad** - Slide off the track of the forehead pad support (Figure 7).
3. **Mask Cushion** - Grasp and pull it off the mask frame (Figure 8).


Assembly of the Mask

1. **Forehead pad** - Slide on the track of the forehead pad support.

▶ Figure 7

2. **Mask Cushion** - Use your thumbs to press the cushion past the retaining tabs on the elbow and press into place (Figure 9).


▶ Figure 8

3. **Headgear** - Thread the top straps through the forehead pad support, and the bottom straps through the headgear clips. Push the headgear clips onto the mask frame.

Tubing Quick Release

The elbow is equipped with a quick release swivel for the CPAP tubing. Simply pull on the swivel with the tubing attached and the tubing will release from the entrainment elbow.


▶ Figure 9

Specifications

⚠ Warning: The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance with Anti-Asphyxia Valve Closed to Atmosphere

Drop in Pressure at

50 SLPM 100 SLPM

All sizes 0.6 cm H₂O 1.5 cm H₂O

Inspiratory and Expiratory Resistance with Anti-Asphyxia Valve Open to Atmosphere

Inspiratory Resistance 0.6 cm H₂O

Expiratory Resistance 0.7 cm H₂O

Deadspace Volume

P 224 ml

S 231 ml

M 242 ml

L 253 ml

Sound Levels

A-weighted Sound Power Level

28 dBA

A-weighted Sound Pressure Level at 1 m

20 dBA

Disposal

Dispose of in accordance with local regulations.

Storage Conditions

Temperature: -4° to 140° F (-20° to +60° C)

Relative Humidity: 15% to 95% non-condensing

Limited Warranty

Respironics, Inc. warrants that its mask systems (including mask frame and cushion) (the "Product") shall be free from defects of workmanship and materials for a period of ninety (90) days from the date of purchase (the "Warranty Period").

If the Product fails under normal conditions of use during the Warranty Period and the Product is returned to Respironics within the Warranty Period, Respironics will replace the Product. This warranty is nontransferable and only applies to the original owner of the Product. The foregoing replacement remedy will be the sole remedy for breach of the foregoing warranty.

This warranty does not cover damage caused by accident, misuse, abuse, negligence, alteration, failure to use or maintain the Product under conditions of normal use and in accordance with the terms of the product literature, and other defects not related to materials or workmanship.

This warranty does not apply to any Product that may have been repaired or altered by anyone other than Respironics. Respironics disclaims all liability for economic loss, loss of profits, overhead, or indirect, consequential, special or incidental damages which may be claimed to arise from any sale or use of the Product. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER EXPRESS WARRANTIES. IN ADDITION, ANY IMPLIED WARRANTY, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IS SPECIFICALLY DISCLAIMED. SOME JURISDICTIONS DO NOT ALLOW DISCLAIMERS OF IMPLIED WARRANTIES, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS UNDER LAWS APPLICABLE IN YOUR SPECIFIC JURISDICTION.

To exercise your rights under this limited warranty, contact your local authorized Respironics, Inc. dealer or Respironics, Inc. at 1001 Murry Ridge Lane, Murrysville, Pennsylvania 15668, USA, or Respironics Deutschland GmbH & Co KG, Gewerbestraße 17, 82211 Herrsching, Germany.

Limited Warranty for Australia and New Zealand only

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above:

Respironics, Inc., a Philips Healthcare company warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 30 days from the date of purchase from an authorised Respironics dealer. 3. If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Respironics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respironics dealer and collecting the product from the authorised Philips Respironics dealer after repair or replacement, at its own cost. Philips Respironics is responsible only for the freight cost of transporting the product between the authorised Philips Respironics dealer and Respironics. Respironics reserves the right to charge an evaluation and postage fee for any returned Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Respironics herein is not transferable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorised Respironics dealer. 7. To the extent permitted by law, where the Buyer has the benefit

of an implied guarantee under the Australian Consumer Law but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Respironics' liability shall be limited at the option of Respironics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty contact your local authorised Philips Respironics dealer. A list of all authorised dealers is available at the following link: http://www.healthcare.philips.com/au_env/homehealth/distributors_index.wpd.

Alternatively, you can make a claim under this warranty by contacting Respironics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde, NSW 2113, Australia. Tel: 1300 766 488, Email: prcontact@philips.com. 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd,
65 Epping Road, North Ryde, NSW 2113
Australia

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©2015 Koninklijke Philips N.V. All rights reserved.